

# An Efficient Technique to Detect Sperm Reactive Oxygen Species: The CellRox Deep Red<sup>®</sup> Fluorescent Probe

Máira Bianchi Rodrigues Alves<sup>1</sup>, André Furugen Cesar de Andrade<sup>2</sup>, Rubens Paes de Arruda<sup>3</sup>, Leonardo Batissaco<sup>1</sup>, Shirley Andrea Florez-Rodriguez<sup>1</sup>, Renata Lançoni<sup>3</sup>, Bruna Marcelle Martins de Oliveira<sup>1</sup>, Mariana Andrade Torres<sup>2</sup>, Gisele Mouro Ravagnani<sup>2</sup>, Tamie Guibu de Almeida<sup>1</sup>, Vinicius Silva Vellone<sup>1</sup> and Eneiva Carla Carvalho Celeghini<sup>1\*</sup>

<sup>1</sup>Laboratory of Research in Pathology of Reproduction, Department of Animal Reproduction, School of Veterinary Medicine and Animal Science, University of Sao Paulo (USP), Pirassununga, SP, Brazil

<sup>2</sup>Laboratory of Andrology and Embryo Technology, Department of Animal Reproduction, School of Veterinary Medicine and Animal Science, University of Sao Paulo (USP), Pirassununga, SP, Brazil

<sup>3</sup>Laboratory of Semen Biotechnology and Andrology, Department of Animal Reproduction, School of Veterinary Medicine and Animal Science, University of Sao Paulo (USP), Pirassununga, SP, Brazil

## Abstract

Heat stress and testicular traumas increase sperm Reactive Oxygen Species (ROS) resulting in Oxidative Stress (OS) and injury of sperm. The fluorescent probe CellROX Deep Red<sup>®</sup> is not already used to detect ROS in spermatozoa. On this way, this study aims to evaluate its efficiency in detecting ROS in ram sperm. For that, it was used ejaculates from rams performing the *in vitro* induction of sperm OS in T0, T50 and T100. T0 was semen sample that was not submitted to OS induction, T50 was 50% induced to OS induction and T100 was entire submitted to OS induction. The production of ROS was detected by CellROX Deep Red<sup>®</sup>. Data polynomial regression analysis was performed and the result of determination coefficient was 0.728. Besides, sixteen rams were submitted to scrotal insulation during 72 hours performing the *in vivo* induction of OS sperm. Analyses of sperm motility, morphology, DNA fragmentation and ROS (detected by CellROX Deep Red<sup>®</sup>) were done two times before the scrotal insulation and two times afterwards. Data analysis of variance was performed from the periods (before x after) and it is observed that after scrotal insulation the total and progressive motility decrease ( $p < 0.05$ ) and total and major sperm defects and sperm DNA fragmentation increase ( $p < 0.05$ ). Moreover, after insulation, it was observed increase ( $p < 0.05$ ) in sperm showing moderate and intense OS. Thus, it is possible to conclude that CellROX<sup>®</sup> fluorescent probe is able to detect ROS production in ram sperm by *in vitro* and *in vivo* induction being a new technique to detect sperm OS.

**Keywords:** Heat stress; Reactive oxygen species; Ram; Male; Infertility; Oxidative stress induction

## Introduction

The global environment temperature is growing and it is already established that it prejudices the male reproduction function. High temperature leads increase of reactive oxygen species (ROS) and, consequently, oxidative stress [1]. Sperm cells are highly susceptible to oxidative stress, which results in fewer ROS neutralizing enzymes in their characteristic small cytoplasm. This leads to a diminished cellular protection. In addition, sperm contain large amounts of polyunsaturated fatty acids, proteins and DNA, which are substrates for ROS [2,3]. Besides, sperm mitochondria produces the most ROS [4], which trigger lipid peroxidation of polyunsaturated fatty acids in the plasma membrane [3,5]. Lipid peroxidation, in turn, reduces plasma membrane fluidity and, as a result, it prejudices sperm function, culminating in decreased sperm motility [6-8] and increased rates of sperm DNA fragmentation [9], cellular apoptosis [2,8,10,11] and morphologically abnormal sperm [2,6].

The detection of ROS in sperm samples is critical for diagnosing infertility. Chemiluminescence is the most commonly used technique for this purpose, but it does not differentiate between ROS produced in sperm and ROS produced in leucocytes. Additionally, results with chemiluminescence have been inconsistent [12]. A widely used technique for measuring oxidative stress, on the other hand, is measuring malondialdehyde with the thiobarbituric acid reaction technique (TBARS) in spectrophotometry [6,13]. However, this technique may produce overestimated values because other substances may also react to TBARS and result in products with an absorbance similar to that of malondialdehyde [14]. Moreover, pH and temperature

alterations can induce peroxidation and produce a reading with higher malondialdehyde values [14].

Various fluorescent probes are available for oxidative stress analysis, such as C11-BODIPY [15], MitoSOX Red<sup>®</sup> [5], and H<sub>2</sub>DCF [16]. C11-BODIPY has affinity for the sperm plasma membrane. It acts as a polyunsaturated fatty acid and, consequently, it is peroxidized by ROS. Thus, in this condition, the fluorescence emitted by C11-BODIPY is altered from red to green [17,18]. Therefore, peroxidation values obtained with this technique are indirect because they express the cell's susceptibility to oxidative stress. On the other hand, Mito SOX Red<sup>®</sup> and H<sub>2</sub>DCF detect the ROS compounds. They are able to detect the superoxide anion and hydrogen peroxide, respectively [16,19]. None of the fluorescent probes used in routine detect hydroxyl radical. On this way, researches of alternative methods for analyzing oxidative stress in sperm are necessary. They must show methods that are easy to apply and interpret.

CellROX Deep Red Reagent<sup>®</sup> probe (CAT 10422 – Molecular Probes) can be analyzed by spectrophotometry, fluorescence microscopy, or flow

**\*Corresponding author:** Eneiva Carla Carvalho Celeghini, Animal Reproduction Department, FMVZ, USP, Av. Duque de Caxias Norte, 225. Postal Code 13635-900, Pirassununga/SP, Brazil, Tel: +55(19) 3565-4242, E-mail: [celeghin@usp.br](mailto:celeghin@usp.br)

Received February 28, 2015; Accepted March 25, 2015; Published April 02, 2015

**Citation:** Alves MBR, de Andrade AFC, de Arruda RP, Batissaco L, Florez-Rodriguez SA, et al. (2015) An Efficient Technique to Detect Sperm Reactive Oxygen Species: The CellRox Deep Red<sup>®</sup> Fluorescent Probe. Biochem Physiol 4: 157. doi: 10.4172/2168-9652.1000157

**Copyright:** © 2015 Alves MBR, et al. This is an open-access article distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited.

cytometry. In addition, it is used in samples preserved in formaldehyde [20] and in combination with other probes [21]. CellROX can be used to measure cellular ROS in the cytoplasm because it detects hydroxyl better than superoxide [22]. It is oxidized during cellular oxidative stress [23] and shows maximum excitation/emission at 640/655 nm. In the absence of ROS, it remains in its reduced state and emits no fluorescence. CellROX is not commonly used in cell biology, and to the best of our knowledge, there are no published studies reporting its use in sperm.

On this way, this study evaluated the efficiency of the use of the CellROX Deep Red Reagent fluorescent probe to detect ROS production in ram sperm. The use of the probe was evaluated in two steps. The first step was performed in experiment 1 that evaluated the efficiency of CellROX to detect sperm oxidative stress induced *in vitro*. In the second experiment, it was evaluated the efficiency of CellROX to detect sperm oxidative stress induced *in vivo*. Both of experiments were done to evaluate the efficiency of CellROX Deep Red Reagent fluorescent probe to detect ROS in ram sperm.

## Materials and Methods

### Experiment 1: Evaluation of CellROX Deep Red Reagent fluorescent probe efficiency in detect sperm

**Oxidative stress induced *in vitro*:** The experiment was conducted with three semen treatments: T0 (semen sample that was not submitted to oxidative stress induction), T50 (semen sample that was 50% induced to oxidative stress and 50% not induced) and T100 (semen sample that was entire submit to oxidative stress induction); and two variables: sperm under mild or no oxidative stress and sperm under moderate or intense oxidative stress.

**Semen collection and preparation:** Three White Dorper rams (*Ovis aries*) with an average age of  $15.3 \pm 2.3$  months and  $66.3 \pm 7.5$  kg were used. The rams were housed in paddock. Hay and concentrate were provided to attend the NRC (1998) [23]. The analyses were conducted on four replicates of semen samples from each of rams ( $n=12$ ). Each of the four replicates was collected as two ejaculates from the same animal harvested within a 90-minute interval; therefore, a total of 24 ejaculates were collected. Samples were collected using an artificial vagina, and only ejaculates showing subjective motility  $\geq 70\%$  and  $\leq 10\%$  total defects were used in the study. The experiment agree with Ethical Principles in Animal Research adopted by "Ethic Committee in the use of animals" of the School of Veterinary Medicine and Animal Science of University of São Paulo, protocol number 2467/2012.

The first ejaculate of each replicate was diluted in TALP media [24] for the preparation of two samples: control samples without oxidative stress induction (T0) and samples with induced oxidative stress (T100). Oxidative stress was induced by adding 50  $\mu$ l of iron sulfate (4 mM) and 50  $\mu$ l of sodium ascorbate (20 mM) to 100  $\mu$ l of semen diluted in TALP and incubated at 37 °C for 90 min. The second ejaculate, collected after the T100 sample incubation, was used for the preparation of the third sample, which contained 50% semen from the second ejaculate diluted in TALP (without oxidative stress induction) and 50% semen from the T100 sample (after oxidative stress induction).

**Oxidative stress analysis:** The CellROX Deep Red Reagent fluorescent probe (2.5 mM; Life Technologies, New York, USA) was diluted in dimethyl sulfoxide (DMSO; 472301, Sigma-Aldrich, St. Louis, USA) for a final concentration of 1 mM (working solution) and

stored at -20 °C in the dark. During use, this working solution was kept in the dark at 37 °C.

Two hundred microliters ( $25 \times 10^6$  sperm/mL) from each semen sample (T0, T50, and T100) were added to 0.5  $\mu$ L of CellROX (1 mM) and 2  $\mu$ L of Hoescht 33342 (in dPBS; 2.5 mg/mL; Life Technologies, New York, USA) and incubated at 37°C for 30 min. After incubation, each sample was centrifuged for 5 minutes at 2,000 g, the supernatant was removed, and the pellet was resuspended in 200  $\mu$ L of TALP.


An aliquot of 4  $\mu$ L of the solution stained with CellROX was placed between a slide and coverslip and read under epifluorescence microscopy (Nikon, model 80i) at 1,000x magnification using a triple filter (D/F/R, C58420) featuring the UV-2E/C (340-380 nm excitation and 435-485 emission), B-2E/C (465-495 excitation and 515-555 emission), and G-2E/C (excitation 540-525 and 605-655 emission) sets.

In each sample was counted 200 cells by the same technique in all analyses. The cells were classified in three categories, as follows: sperm under mild or no oxidative stress (unstained midpiece), sperm under moderate oxidative stress (midpiece stained pale red), and sperm under intense oxidative stress (midpiece stained strong red), as showed in the Figure 1.

**Statistical analysis:** We conducted a polynomial regression analysis using the Statview software (Stat view 1998, SAS Institute Inc., Cary, NC, USA). The experimental model was:  $Y = a + bX$ , wherein  $Y$  is the estimative of the oxidative stress as a function of the treatment,  $a$  is the polynomial regression coefficient corresponding to the value of  $Y$  when  $X$  is 0,  $b$  is the regression coefficient for the percent of  $X$  on the response  $Y$ , and  $X$  is the treatment.

### Experiment 2: Oxidative stress induction by testicular degeneration in ram sperm

**Animals:** Semen samples from 16 White Dorper rams with an average age of  $17.62 \pm 2.6$  months and  $67.12 \pm 7.2$  kg were analyzed. The rams were housed in paddock. Hay and concentrate were provided to attend the NRC (1998). Prior to the beginning of the study, sperm motility was  $\geq 60\%$  and total defects were  $\leq 20\%$ . The animals were submitted to scrotal insulation for 72 hours through the use of heating bags to induce testicular degeneration and loss of testicular thermoregulation. We performed four semen collections with 15-day intervals in between; two were performed before the scrotal insulation and two afterwards. Semen samples were collected using an artificial vagina. The experiment agree with Ethical Principles in Animal



**Figure 1:** Epifluorescence photomicrography of ram sperm cells stained with CellROX Deep Red® and Hoescht 33342 fluorescent probes. **A:** sperm under mild or no oxidative stress (middle piece no stained). **B:** sperm under moderate oxidative stress (middle piece stained in pale red). **C:** sperm under intense oxidative stress (middle piece stained in intense red).

Research adopted by "Ethic Committee in the use of animals" of the School of Veterinary Medicine and Animal Science of University of São Paulo, protocol number 2467/2012.

**Semen preparation and evaluation of testicular degeneration:** For examination of the testicular degeneration induction was performed evaluation of total and progressive motility, sperm defects and sperm DNA fragmentation. The total and progressive motility was assessed by computer analyses system. The software was the Sperm Class Analyzer (SCA, Microoptics, Spain) and the set up was adjusted by ram sperm. It was chosen five fields for the analyses and it was used the Makler chamber to perform the evaluation. The sperm defects were evaluated by differential interference contrast microscopy (DIC). It was counted 200 cells that were classified in major and minor sperm defects. DNA fragmentation was determined by commercial kit Halomax (Halotech, Spain) and it was determined the fragmentation in 1,000 cells. All evaluation was done by competent techniques.

**Semen preparation and evaluation of oxidative stress:** Two hundred microliters of semen samples diluted in TALP media [25] ( $25 \times 10^6$  sperm/mL) were added to 0.5  $\mu$ L of CellROX (1 mM in DMSO) and 2  $\mu$ L of Hoescht 33342 (2.5 mg/mL; in dPBS), and incubated for 30 minutes at 37°C. After incubation, the solution was centrifuged for 5 minutes at 2,000 g, the supernatant was removed, and the pellet resuspended in 200  $\mu$ L of TALP. An aliquot of 5  $\mu$ L was placed between a slide and coverslip and read under epifluorescence microscopy (Nikon, model 80i) at 1,000x magnification using a triple filter (D/F/R, C58420) featuring the UV-2E/C (340-380 nm excitation and 435-485 emission), B-2E/C (465-495 excitation and 515-555 emission), and G-2E/C (540-525 excitation and 605-655 emission) sets.

It was counted 200 cells and they were classified as sperm under mild or no oxidative stress (unstained midpiece), sperm under moderate oxidative stress (midpiece stained pale red), and sperm under intense oxidative stress (midpiece stained strong red) (Figure 1).


**Statistical analysis:** Nine variables were analyzed in this experiment (total motility, progressive motility, total sperm defects, minor sperm defects, major sperm defects, DNA fragmentation, sperm under mild or no oxidative stress, sperm under moderate oxidative stress, and sperm under intense oxidative stress) in two moments (before insulation/control period, and after insulation/degenerated period).

An analysis of variance was conducted (ANOVA). Differences between treatments were considered significant at  $p < 0.05$  (Stat view 1998, SAS Institute Inc., Cary, NC, USA).

## Results

In the first experiment, the percentage of cells with moderate or intense oxidative stress was compared between treatments T0, T50 and T100. T0, T50 and T100 showed, respectively,  $0.25 \pm 0.11\%$ ,  $16.37 \pm 1.75\%$  and  $30.83 \pm 3.57\%$  of moderate or intense oxidative stress. Results from the polynomial regression analysis ( $R^2=0.728$ ) are shown in Figure 2.

In the second experiment, the testicular degeneration and oxidative stress were induced by scrotal insulation and it caused a decrease in total ( $p < 0.0001$ ) and progressive motility ( $p < 0.0001$ ), and an increase in total ( $p < 0.0001$ ) and major sperm ( $p < 0.0001$ ) defects and DNA fragmentation ( $p=0.001$ ) (Table 1). Besides, it caused a reduction in sperm with mild or no oxidative stress ( $p=0.02$ ), and an increase in sperm showing moderate ( $p=0.03$ ) and intense ( $p=0.03$ ) oxidative stress (Table 2).



**Figure 2:** Polynomial regression analysis of oxidative stress detected in ram sperm with the Cell ROX Deep Red® fluorescent probe submitted to the T0 (without oxidative stress induction), T50 (50% of the sample without and 50% of the sample with oxidative stress induction), and T100 (with oxidative stress induction) treatments. It was considered spermatozoa with mild (cells with middle piece no stained), moderate (cells with middle piece stained in pale red) and intense oxidative stress (cells with middle piece stained in intense red).

Semen evaluation (%)	Control	Testicular Degeneration	P value
Total motility	87.25 ± 1.31	45.63 ± 5.01	<0.0001
Progressive motility	59.59 ± 2.10	30.64 ± 3.89	<0.0001
Total defects	19.24 ± 2.36	64.46 ± 3.40	<0.0001
Major defects	9.69 ± 1.37	55.46 ± 3.73	<0.0001
Minor defects	9.55 ± 1.73	9.00 ± 2.19	0.8457
DNA fragmentation	1.10 ± 0.28	8.50 ± 2.20	0.0017

**Table 1:** Average ± standard error of total motility (%), progressive motility (%), total (%), major (%) and minor (%) sperm defects and DNA fragmentation (%) before (control) and after testicular degeneration in rams.

Degree of Oxidative Stress (%)	Control	Testicular Degeneration	P value
Absent or mild oxidative stress	84.10 ± 2.94	73.75 ± 3.56	0.0287
Moderate oxidative stress	15.68 ± 2.92	25.43 ± 3.41	0.0339
Intense oxidative stress	0.20 ± 0.13	0.81 ± 0.24	0.0343

**Table 2:** Average ± standard error of oxidative stress in sperm (%) before (control) and after testicular degeneration in rams.

## Discussion

The use of the CellROX Deep Red® probe in the analysis of sperm oxidative stress using ovine species as the experimental model was evaluated and it was successfully validated for the detection of oxidative stress induction. This fluorescent probe has not been previously used in sperm cells.

Ferrous sulfate ( $\text{FeSO}_4$ ) and sodium ascorbate were used to induce oxidative stress *in vitro* because a preliminary test using menadione, cited by the commercial guide as an inducer of oxidative stress in bovine lung endothelial cells [26], did not effectively induce oxidative stress in ram sperm. In boar sperm, the induction of oxidative stress with  $\text{FeSO}_4$  and Na ascorbate was also more efficient than with menadione and hydrogen peroxide [17]. After 30 minutes' incubation with  $\text{FeSO}_4$  and Na ascorbate, an analysis with BODIPY showed  $92.9 \pm 5.9\%$  oxidized sperm, which was significantly different from the non-induced sample ( $p < 0.05$ ). After 30 minutes' incubation with menadione and hydrogen peroxide, the peroxidation rates were  $0.3 \pm 0.1\%$  and  $1.1 \pm 0.3\%$ , respectively. These results were not significantly different from the non-induced sample [17].

In the present study, the rates of CellROX® positive sperm observed

in the T100 group were lower than those reported by Guthrie and Welch (2007) [10]. These differences may be explained by the differences between the CellROX<sup>®</sup> and BODIPY methods. CellROX<sup>®</sup> assesses the amount of hydroxyl and peroxide present in the sample. BODIPY, used by Guthrie and Welch (2007) [10], assesses the lipid peroxidation caused by any ROS [18]. We know that antioxidant agents in the seminal plasma [27,28] neutralize the ROS molecules present there [3]. On this way, the BODIPY assesses more ROS than CellROX<sup>®</sup>. Nevertheless, regardless of the neutralizing action of the seminal plasma, the amounts of ROS (including non-neutralized) in T50 and T100 were sufficient to produce different results between groups.

The rams in experiment 1 were healthy and did not present signs of testicular degeneration, showing  $\geq 70\%$  sperm motility, and  $\leq 10\%$  sperm defects. Thus, the ROS rates were expected to be lower in T0 and T50 than in T100. The same result was expected for the two retrievals conducted prior to the use of heating bags.

The loss of testicular thermoregulation was induced by applying heating bags to the animals' testes. The confirmation of thermoregulation deleterious was done evaluating the motility, sperm defects and DNA fragmentation. Loss of testicular thermoregulation leads to increase cellular metabolism, hypoxia, and a subsequent increase in ROS production, which in turn leads to oxidative stress [8,11,29,30]. Mice subjected to thermal stress at 42°C for 30 minutes showed an increased expression of heme oxygenase-1 after 24 and 48 hours of stress exposure, and increased glutathione peroxidase 1 and alpha glutathione-S-transferase after 3 hours [11]. A significant increase in 4-HNE protein expression was observed in men diagnosed with different degrees of varicocele compared to those in the control group, while a positive significant correlation ( $p=0.003$ ;  $r=0.49$ ) between the increased temperature and 4-HNE protein expression was observed in the groups with varicocele [31]. Therefore, after the insulation period, we expected to see an increase in sperm cell ROS production, as we previously saw for the last two stages.

On this way, the results showed that the CellROX Deep Red<sup>®</sup> fluorescent probe is efficient in detecting the production of sperm ROS under oxidative stress induced *in vitro* and *in vivo*. The detection of ROS in ram sperm can be important in the diagnosis of reproductive disorders.

#### Acknowledgements

The authors thank Dr. Marcílio Nichi and Dr. João Diego de Agostini for assistance in the performance of the oxidative stress induction. We also thank Mr. José Maria Bernardes, Márcio De Carli and João Carlos Campos for assistance to the animals. This work was supported by the Fundação de Amparo à Pesquisa do Estado de São Paulo (FAPESP process numbers 2011/16744-3 and 2012/00040-0).

#### References

- Hansen PJ (2009) Effects of heat stress on mammalian reproduction. *Philos Trans R Soc Lond B Biol Sci* 364: 3341-3350.
- Aitken RJ, Jones KT, Robertson SA (2012) Reactive oxygen species and sperm function—in sickness and in health. *J Androl* 33: 1096-1106.
- Walczak-Jedrzejowska R, Wolski JK, Slowikowska-Hilczner J (2013) The role of oxidative stress and antioxidants in male fertility. *Cent European J Urol* 66: 60-67.
- Koppers AJ, Garg ML, Aitken RJ (2010) Stimulation of mitochondrial reactive oxygen species production by unesterified, unsaturated fatty acids in defective human spermatozoa. *Free Radic Biol Med* 48: 112-119.
- Aitken RJ, Smith TB, Lord T, Kuczera L, Koppers AJ, et al. (2013) On methods for the detection of reactive oxygen species generation by human spermatozoa: analysis of the cellular responses to catechol oestrogen, lipid aldehyde, menadione and arachidonic acid. *Andrology* 1: 192-205.
- Aitken RJ, Harkiss D, Buckingham D (1993) Relationship between iron-catalysed lipid peroxidation potential and human sperm function. *J Reprod Fertil* 98: 257-265.
- Aitken RJ, De Luliis GN, Gibb Z, Baker MA (2012) The Simmet lecture: new horizons on an old landscape—oxidative stress, DNA damage and apoptosis in the male germ line. *Reprod Domest Anim* 47 Suppl 4: 7-14.
- Pérez-Crespo M, Pintado B, Gutiérrez-Adán A (2008) Scrotal heat stress effects on sperm viability, sperm DNA integrity, and the offspring sex ratio in mice. *Mol Reprod Dev* 75: 40-47.
- Paul C, Murray AA, Spears N, Saunders PT (2008) A single, mild, transient scrotal heat stress causes DNA damage, subfertility and impairs formation of blastocysts in mice. *Reproduction* 136: 73-84.
- Guthrie HD, Welch GR (2012) Effects of reactive oxygen species on sperm function. *Theriogenology* 78: 1700-1708.
- Paul C, Teng S, Saunders PT (2009) A single, mild, transient scrotal heat stress causes hypoxia and oxidative stress in mouse testes, which induces germ cell death. *Biol Reprod* 80: 913-919.
- Aitken RJ, West KM (1990) Analysis of the relationship between reactive oxygen species production and leucocyte infiltration in fractions of human semen separated on Percoll gradients. *Int J Androl* 13: 433-451.
- Nichi M, Bols PE, Züge RM, Barnabe VH, Goovaerts IG, et al. (2006) Seasonal variation in semen quality in *Bos indicus* and *Bos taurus* bulls raised under tropical conditions. *Theriogenology* 66: 822-828.
- Janero DR (1990) Malondialdehyde and thiobarbituric acid-reactivity as diagnostic indices of lipid peroxidation and peroxidative tissue injury. *Free Radic Biol Med* 9: 515-540.
- Collodel G, Federico MG, Geminiani M, Martini S, Bonechi C, et al. (2011) Effect of trans-resveratrol on induced oxidative stress in human sperm and in rat germinal cells. *Reprod Toxicol* 31: 239-246.
- Guthrie HD, Welch GR (2006) Determination of intracellular reactive oxygen species and high mitochondrial membrane potential in Percoll-treated viable boar sperm using fluorescence-activated flow cytometry. *J Anim Sci* 84: 2089-2100.
- Guthrie HD, Welch GR (2007) Use of fluorescence-activated flow cytometry to determine membrane lipid peroxidation during hypothermic liquid storage and freeze-thawing of viable boar sperm loaded with 4-(difluoro-5-(4-phenyl-3-butadienyl)-4-bora-3a,4a-diaza-s-indacene-3-undecanoic acid. *J Anim Sci* 85: 1402-1411.
- Pap EH, Drummen GP, Winter VJ, Kooij TW, Rijken P, et al. (1999) Ratio-fluorescence microscopy of lipid oxidation in living cells using C11-BODIPY(581/591). *FEBS Lett* 453: 278-282.
- Koppers AJ, De Luliis GN, Finnie JM, McLaughlin EA, Aitken RJ (2008) Significance of mitochondrial reactive oxygen species in the generation of oxidative stress in spermatozoa. *J Clin Endocrinol Metab* 93: 3199-3207.
- Grinberg YY, Dibbern ME, Levasseur VA, Kraig RP (2013) Insulin-like growth factor-1 abrogates microglial oxidative stress and TNF- $\alpha$  responses to spreading depression. *J Neurochem* 126: 662-672.
- Ahn HY, Fairfull-Smith KE, Morrow BJ, Lussini V, Kim B, et al. (2012) Two-photon fluorescence microscopy imaging of cellular oxidative stress using profluorescent nitroxides. *J Am Chem Soc* 134: 4721-4730.
- Grinberg Y (2013) ROS production measurement using con-focal microscopy?
- Grinberg YY, van Drongelen W, Kraig RP (2012) Insulin-like growth factor-1 lowers spreading depression susceptibility and reduces oxidative stress. *J Neurochem* 122: 221-229.
- National Research Council (1984) Nutrient Requirement of Beef Cattle. National Academies Press, Washington.
- Bavister D, Lorraine M (1983) Development of preimplantation embryos of the golden hamster in a defined culture medium. *Biol Reprod* 28: 235-247.
- The Molecular Probes<sup>®</sup> Handbook - A Guide to Fluorescent Probes and Labeling Technologies. (11thedn).
- Li TK (1975) The glutathione and thiol content of mammalian spermatozoa and seminal plasma. *Biol Reprod* 12: 641-646.
- Potts RJ, Notarianni LJ, Jefferies TM (2000) Seminal plasma reduces exogenous

- oxidative damage to human sperm, determined by the measurement of DNA strand breaks and lipid peroxidation. *Mutat Res* 447: 249–256.
29. Hansen PJ (2009) Effects of heat stress on mammalian reproduction. *Philos Trans R Soc Lond B Biol Sci* 364: 3341–3350.
30. Setchell BP (1998) The Parkes Lecture. Heat and the testis. *J Reprod Fertil* 114: 179–194.
31. Shiraishi K, Takihara H, Matsuyama H (2010) Elevated scrotal temperature, but not varicocele grade, reflects testicular oxidative stress-mediated apoptosis. *World J Urol* 28: 359–364.