4th International Conference and Expo on

## **Novel Physiotherapies**

August 21-22, 2017 | Birmingham, UK

## Treatment of supraspinatus tendinopathy with Ergon® IASTM technique and neuromuscular control exercises: A case study

Konstantinos Fousekis, Konstantinos Mylonas and Pavlos Angelopoulos Technological Educational Institute of Western Greece, Greece

**Study Background:** Supraspinatus tendinopathy is an important cause of pain and dysfunction in the adult shoulder. Traditional treatment for this type of injury includes traditional forms of treatment such as massage, stretching, electrotherapy and functional exercises. The aim of this case study was to evaluate the effectiveness of Ergon\* IASTM Technique and shoulder neuromuscular control exercises in the treatment of supraspinatus tendinopathy.

**Methods:** A 50-year-old patient clinically diagnosed with supraspinatus tendinosis presented with significant a) pain on palpation b) pain during passive and active internal rotation and c) decrease in shoulder passive internal rotation. His treatment plan included 8 treatment sessions involving the application a) of Ergon® IASTM Technique over specific shoulder points and b) of targeted neuromuscular control exercises of the shoulder. Pain produced during passive internal rotation was evaluated with a visual analogue scale (VAS). The range of motion (ROM) for the internal rotation was measured with a goniometer. The patient was evaluated before, and at the 4th and 8th treatment session.

**Results:** The patient experienced a significant decrease in pain and an increase in shoulder ROM regarding internal rotation by both the 4th and 8th treatment session (figure1). More specifically, pain, as measured by VAS scale, was decreased from 8 and 7, respectively on the passive and active internal rotation of the shoulder, to 6 and 4 by the end of the 4th week and to 3 and 2 after the 8th treatment. Internal rotation ROM in the painful shoulder at 90° of abduction progressed from 60° at the baseline to 730 and 780 after 4th and 8th treatment, respectively.

**Conclusions:** This case study provides some evidence that Ergon® IASTM Technique in association with shoulder neuromuscular control exercises is an effective technique in the rehabilitation of the patients with supraspinatus tendinopathy.

## **Biography**

Konstantinos Fousekis is an Associate Professor in Sports Physiotherapy at the Department of Physical Therapy, Technological Educational Institute of Western Greece. He is a Physiotherapist specializing in soft tissue mobilization techniques (IASTM). He has years of experience in treating musculoskeletal and sports injuries and is a Professional Physical Therapist for several professional soccer teams. His research interests deals with the assessment and rehabilitation of sports and musculoskeletal injuries using IASTM techniques. In cooperation with Konstantinos Mylonas, he created the ERGON® IASTM Technique as a basic treatment of painful and non-musculoskeletal disorders.

konfousekis@gmail.com

TIAN T			
	Ot	OC	0
Τ.4	υı	CO	۰